

# Erba Chem 5x<sup>®</sup>

- Extensive Biochemistry test menu
- Elaborate Q.C. Program and Memory
- Inbuilt Printer
- Graphic Display
- 1000 test Result storage



**Clinical Chemistry Analyzer**



# Erba Chem 5<sup>®</sup>X


Clinical Chemistry Analyzer

## TECHNICAL SPECIFICATIONS

- 150 open user-defined test programs
- 9 Assay modes : 1 point Linear, 2 point Linear, 1 point Non Linear, 2 point Non-linear, Rate A Linear, Rate A Non-Linear, 1 point Sample Blank Linear, 1-point Sample Blank Non-Linear and Absorbance
- Mono & Bichromatic measurement
- Photometric range from 0 to 3.0 OD
- Photometric resolution up to 0.0001
- Automatic Zero Setting
- Static Photometer (340-630 nm)
- Wavelength 340, 405, 505, 546, 578, 630 nm
- Multi-Point calibration with graphical display
- Levy-Jennins graph for quality control
- 1000 Test Patient Results storage
- Peristaltic Pump Aspiration System with calibration facility
- Dual Cuvette System : Flow cell volume-33  $\mu$ l  
: 10 mm external glass / Polystyrene
- Aspiration Volume from 350 to 999  $\mu$ l
- Temperature Control at 25°, 30° & 37° C
- Repeat measurement by Read key
- Collated report date wise, ID wise (Alphanumeric)
- Inbuilt Thermal Printer
- Direct Printer Connectivity through USB Port
- USB Connectivity to PC
- Connectivity to external Incubator ( Optional with 16 position )

Specifications are subject to change without prior notification.

**Distributor Stamp**



**ERBA Diagnostics Mannheim GmbH**

Mallastrasse 69-73 68219 Mannheim, Germany  
Telephone : (+49) 621 8799770 Fax : (+49) 621 8799688  
Email : sales@erbamannheim.com Website : www.erbamannheim.com

**TRANSASIA BIO - MEDICALS LTD.**

Transasia House, 8 Chandivali Studio Road, Andheri (E), Mumbai - 400 072  
Tel.: (022) 4030 9000, Fax : (022) 2857 3030 Email : transasia@transasia.co.in  
Website : [www.transasia.co.in](http://www.transasia.co.in) Twitter : [https://twitter.com/transasia\\_1](https://twitter.com/transasia_1)